




# UK Standards for Microbiology Investigations

## Coagulase test



"NICE has renewed accreditation of the process used by **Public Health England (PHE)** to produce **UK Standards for Microbiology Investigations**. The renewed accreditation is valid until **30 June 2021** and applies to guidance produced using the processes described in **UK standards for microbiology investigations (UKSMIs) Development process, S9365', 2016**. The original accreditation term began in **July 2011**."

## Acknowledgments

UK Standards for Microbiology Investigations (UK SMIs) are developed under the auspices of Public Health England (PHE) working in partnership with the National Health Service (NHS), Public Health Wales and with the professional organisations whose logos are displayed below and listed on the website <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. UK SMIs are developed, reviewed and revised by various working groups which are overseen by a steering committee (see <https://www.gov.uk/government/groups/standards-for-microbiology-investigations-steering-committee>).

The contributions of many individuals in clinical, specialist and reference laboratories who have provided information and comments during the development of this document are acknowledged. We are grateful to the medical editors for editing the medical content.

For further information please contact us at:

Standards Unit  
National Infection Service  
Public Health England  
61 Colindale Avenue  
London NW9 5EQ

E-mail: [standards@phe.gov.uk](mailto:standards@phe.gov.uk)

Website: <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>

PHE publications gateway number: 201703

UK Standards for Microbiology Investigations are produced in association with:



Logos correct at time of publishing.

# Contents

---

Acknowledgments .....	2
Contents .....	3
Amendment table .....	4
UK SMI: scope and purpose .....	5
Scope of document .....	7
Introduction .....	7
Technical information/limitations .....	7
1 Safety considerations .....	10
2 Reagents and equipment .....	10
3 Quality control organisms .....	10
4 Procedure and results .....	11
Appendix: Coagulase test .....	14
References .....	15



"NICE has renewed accreditation of the process used by **Public Health England (PHE)** to produce **UK Standards for Microbiology Investigations**. The renewed accreditation is valid until **30 June 2021** and applies to guidance produced using the processes described in **UK standards for microbiology investigations (UKSMIs) Development process, S9365', 2016**. The original accreditation term began in **July 2011**."

## Amendment table

---

Each UK SMI method has an individual record of amendments. The current amendments are listed on this page. The amendment history is available from [standards@phe.gov.uk](mailto:standards@phe.gov.uk).

New or revised documents should be controlled within the laboratory in accordance with the local quality management system.

Amendment number/date	9/01.08.18
Issue number discarded	5
Insert issue number	6
Anticipated next review date*	01.08.21
<b>Section(s) involved</b>	<b>Amendment</b>
Whole document.	<p>Document updated.</p> <p>Technical limitations updated with subheadings.</p> <p>Information on plasma use and NCTC 6571 added to the technical limitations/information.</p> <p>Link to UK SMI TP 34: Thermonuclease test document added to this document.</p> <p>References updated with grades.</p>

\*Reviews can be extended up to five years subject to resources available.

## UK SMI<sup>#</sup>: scope and purpose

---

### Users of UK SMIs

Primarily, UK SMIs are intended as a general resource for practising professionals operating in the field of laboratory medicine and infection specialties in the UK. UK SMIs also provide clinicians with information about the available test repertoire and the standard of laboratory services they should expect for the investigation of infection in their patients, as well as providing information that aids the electronic ordering of appropriate tests. The documents also provide commissioners of healthcare services with the appropriateness and standard of microbiology investigations they should be seeking as part of the clinical and public health care package for their population.

### Background to UK SMIs

UK SMIs comprise a collection of recommended algorithms and procedures covering all stages of the investigative process in microbiology from the pre-analytical (clinical syndrome) stage to the analytical (laboratory testing) and post analytical (result interpretation and reporting) stages. Syndromic algorithms are supported by more detailed documents containing advice on the investigation of specific diseases and infections. Quality guidance notes describe laboratory processes which underpin quality, for example assay validation.

Standardisation of the diagnostic process through the application of UK SMIs helps to assure the equivalence of investigation strategies in different laboratories across the UK and is essential for public health surveillance, research and development activities.

### Equal partnership working

UK SMIs are developed in equal partnership with PHE, NHS, Royal College of Pathologists and professional societies. The list of participating societies may be found at <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. Inclusion of a logo in an UK SMI indicates participation of the society in equal partnership and support for the objectives and process of preparing UK SMIs. Nominees of professional societies are members of the Steering Committee and working groups which develop UK SMIs. The views of nominees cannot be rigorously representative of the members of their nominating organisations nor the corporate views of their organisations. Nominees act as a conduit for two way reporting and dialogue. Representative views are sought through the consultation process. UK SMIs are developed, reviewed and updated through a wide consultation process.

### Quality assurance

NICE has accredited the process used by the UK SMI working groups to produce UK SMIs. The accreditation is applicable to all guidance produced since October 2009. The process for the development of UK SMIs is certified to ISO 9001:2008. UK SMIs represent a good standard of practice to which all clinical and public health microbiology laboratories in the UK are expected to work. UK SMIs are NICE accredited and represent neither minimum standards of practice nor the highest level

---

<sup>#</sup> Microbiology is used as a generic term to include the two GMC-recognised specialties of Medical Microbiology (which includes Bacteriology, Mycology and Parasitology) and Medical Virology.

of complex laboratory investigation possible. In using UK SMIs, laboratories should take account of local requirements and undertake additional investigations where appropriate. UK SMIs help laboratories to meet accreditation requirements by promoting high quality practices which are auditable. UK SMIs also provide a reference point for method development. The performance of UK SMIs depends on competent staff and appropriate quality reagents and equipment. Laboratories should ensure that all commercial and in-house tests have been validated and shown to be fit for purpose. Laboratories should participate in external quality assessment schemes and undertake relevant internal quality control procedures.

## Patient and public involvement

The UK SMI working groups are committed to patient and public involvement in the development of UK SMIs. By involving the public, health professionals, scientists and voluntary organisations the resulting UK SMI will be robust and meet the needs of the user. An opportunity is given to members of the public to contribute to consultations through our open access website.

## Information governance and equality

PHE is a Caldicott compliant organisation. It seeks to take every possible precaution to prevent unauthorised disclosure of patient details and to ensure that patient-related records are kept under secure conditions. The development of UK SMIs is subject to PHE Equality objectives <https://www.gov.uk/government/organisations/public-health-england/about/equality-and-diversity>.

The UK SMI working groups are committed to achieving the equality objectives by effective consultation with members of the public, partners, stakeholders and specialist interest groups.

## Legal statement

While every care has been taken in the preparation of UK SMIs, PHE and the partner organisations, shall, to the greatest extent possible under any applicable law, exclude liability for all losses, costs, claims, damages or expenses arising out of or connected with the use of an UK SMI or any information contained therein. If alterations are made by an end user to an UK SMI for local use, it must be made clear where in the document the alterations have been made and by whom such alterations have been made and also acknowledged that PHE and the partner organisations shall bear no liability for such alterations. For the further avoidance of doubt, as UK SMIs have been developed for application within the UK, any application outside the UK shall be at the user's risk.

The evidence base and microbial taxonomy for the UK SMI is as complete as possible at the date of issue. Any omissions and new material will be considered at the next review. These standards can only be superseded by revisions of the standard, legislative action, or by NICE accredited guidance.

UK SMIs are Crown copyright which should be acknowledged where appropriate.

## Suggested citation for this document

Public Health England. (2018). Coagulase test. UK Standards for Microbiology Investigations. TP 10 Issue 6. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>

## Scope of document

---

This test detects the ability of bacteria to clot plasma by the action of the enzyme, coagulase. This test is used specifically to differentiate the species within the genus *Staphylococcus*. There are many hypotheses but the precise mechanism of coagulase action is not known<sup>1</sup>.

This UK SMI should be used in conjunction with other UK SMIs.

## Introduction

---

Coagulase is a protein enzyme produced by several microorganisms that enables the conversion of fibrinogen to fibrin. Coagulase binds plasma fibrinogen, causing the organisms to agglutinate or plasma to clot. Coagulase exists in two forms: “bound coagulase” (or clumping factor) which is bound to the cell wall and “free coagulase” which is liberated by the cell wall. Bound coagulase is detected by the slide coagulase test, whereas both free and bound coagulase are detected by the tube coagulase test.

Bound coagulase adsorbs fibrinogen from the plasma and alters it so it precipitates on the staphylococci, causing them to clump resulting in cell agglutination. Free coagulase reacts with a substance in plasma to form a fibrin clot.

## Technical information/limitations

---

### Purity of inoculum used

The colony inoculum used for testing must be pure because a contaminant may produce false results after prolonged incubation.

### Slide coagulase test

This test is unsuitable for isolates that are not easily emulsified<sup>2</sup>.

Autoagglutination may occur.

Use water instead of saline as some staphylococci are salt sensitive, particularly if they have been cultured on salt containing media (such as Mannitol Salt Agar), as lysis or clumping of cells may occur. Do not perform slide coagulase tests from colonies grown from salt containing media.

Over mixing may cause the clots to break down<sup>3</sup>.

*S. schleiferi* and *S. lugdunensis* may give positive results on the slide coagulase test<sup>2,4</sup>.

### Tube coagulase test

For the tube coagulase test, EDTA plasma is superior to citrated plasma because citrate-utilizing organisms such as *Pseudomonas* species, *Serratia marcescens*, *Enterococcus faecalis* and strains of *Streptococcus* will clot citrated plasma<sup>5</sup>.

Longer incubation at 37°C may result in disappearance of the clot. This is due to the production of staphylokinase which can lyse the clot<sup>5</sup>.

Some other species of staphylococci, including *Staphylococcus schleiferi* and *Staphylococcus intermedius* may give positive results in the tube coagulase test but are not common isolates from human infections<sup>2,4</sup>.

The tube coagulase test should not be unduly agitated as this can cause the clot to shrink also giving a false negative result<sup>6</sup>.

### Plasma use

Commercially available lyophilized products (or kits) are available and manufacturers' instructions should be followed when using these products.

### Tube coagulase test direct from blood culture broth

The main advantage of the direct tube coagulase test is its rapidity in presumptive identification of *S. aureus* in blood culture broths and reducing turnaround times and it is also helpful in initiating appropriate antimicrobial treatment<sup>7,8</sup>. It is a valuable adjunct in the routine microbiology laboratory because of its good performance, technical simplicity and low cost.

Care should be taken when using this test directly on presumptive positive coagulase blood culture broths. Several investigators have evaluated the accuracy of two agglutination tests in rapid identification of *S. aureus* from positive blood culture broths and although overall specificity has been excellent, a wide range of sensitivities have been reported with a variety of latex tests<sup>7,9,10</sup>. Studies have confirmed that slide coagulase kits should not be used for the rapid identification of *S. aureus* from blood culture broths as they have been principally designed for isolates grown on solid culture media and lack sensitivity. Manufacturer's instruction should be followed.

The tube coagulase test is the preferred method for the rapid screening test for *S. aureus* in blood culture broths and when reported negative, blood culture isolates should be re-tested using standard laboratory techniques (Gram stain, subculturing of broths and retesting from the culture plates). There have been reports of negative coagulase test results when tested directly on blood cultures, but when repeated from culture isolates, are reported positive<sup>7</sup>. Another limitation that may deter laboratories from its use is the labour intensive, double centrifugation step. Moreover, two recent reports have indicated no loss of sensitivity when the tube coagulase test is performed directly on uncentrifuged blood culture broths<sup>10,11</sup>.

Some strains of Meticillin Resistant *Staphylococcus aureus* may exhibit a negative or weak positive reaction. In addition, rare strains of *S. aureus* are negative in coagulase tests<sup>2</sup>.

Refer to UK SMI [TP 34: Thermonuclease test](#) procedure for information on how to use this test to determine the presence of *S. aureus* in positive blood culture bottles.

### Commercial kits

Commercial kits are available using latex technology. These kits can detect Protein A and/or clumping factor but can also detect various surface antigens, making them more sensitive than the coagulase test but at some expense to specificity due to cross-reaction with Coagulase Negative Staphylococci. In addition, any test including clumping factor may give false positive results with *Staphylococcus lugdunensis* and *Staphylococcus schleifer*<sup>2</sup>.



**NCTC 6571**

The NCTC strain known as the “*Oxford Staphylococcus*” is widely used in clinical diagnostic microbiology laboratories throughout the UK as a reference control strain for antimicrobial susceptibility testing and other phenotypic tests such as coagulase test (regarded as a weak positive control), determination of DNase activity, etc.

Recent studies have shown that this strain produces Pantone-Valentine Leukocidin (PVL), a pore-forming cytotoxin produced by fewer than 5% of *Staphylococcus aureus* strains that causes leukocyte destruction and tissue necrosis. Therefore it is recommended that good practice should be adhered to when this organism is handled<sup>12</sup>.

## 1 Safety considerations<sup>13-30</sup>

---

Refer to current guidance on the safe handling of all organisms and reagents documented in this UK SMI.

Aseptic technique and established precautions against microbiological hazards throughout all procedures must be observed.

All work likely to generate aerosols must be performed in a microbiological safety cabinet.

The above guidance should be supplemented with local COSHH and risk assessments.

Compliance with postal and transport regulations is essential.

## 2 Reagents and equipment

---

Fresh discrete bacterial colonies growing on solid medium (non-salt containing media) including the positive and negative control organisms

Distilled water

Microscope slide

Bacteriological loop (preferably nichrome) or disposable alternative

Disposable Pasteur pipette

Marking pen or wax pencil

Sterile test tubes

Test tube rack

### Test solution

#### Slide coagulase test:

Commercially available plasma treated with ethylene diamine-tetraacetic acid (EDTA).

#### Tube coagulase test:

Commercially available plasma treated with EDTA added suitable for tube coagulase. Use the plasma according to manufacturer's instructions unless an alternative method has been validated.

**Note:** Commercially lyophilized products (or kits) are also available and manufacturer's instructions should be followed.

## 3 Quality control organisms

---

Refer to [TP 1- Example reference strains for UK SMI test procedures](#)

### Positive control

*Staphylococcus aureus*

NCTC 6571 or NCTC 12973

## Negative control

*Staphylococcus haemolyticus* NCTC 11042

**Note:** These strains have been validated by NCTC to give this result.

## 4 Procedure and results

### 4.1 Slide coagulase test<sup>1,31</sup>

- place two drops of distilled water on a clean glass slide. Identify where the test strain (T) and the control (C) will be placed by labelling the slide. An additional slide will be required for the control strains and this should be clearly labelled.
- set up the positive and negative control organisms on the same slide to be tested simultaneously.
- emulsify the test strain to obtain a homogenous thick suspension. False negative reactions will occur if the bacterial suspension is not heavy enough.
- observe for autoagglutination. Strains which autoagglutinate must be tested by an alternative procedure.
- dip a straight wire or loop in the plasma and stir gently with the homogenous suspension. If using a reusable loop, sterilize the loop before proceeding.

**Note:** Plasma is added only to the test strain and the control organisms but not the control (C) as it serves as an autoagglutination control.

- observe for immediate formation of white clumps.

#### Positive result

Visible clumping within 10s.

#### Negative result

No visible clumping within 10s.

**Note:** The positive control species should show clumping only when emulsified in the plasma and the negative control species should not show clumping in either water (saline) or plasma<sup>4</sup>.

### 4.2 Tube coagulase test<sup>1-3,6,31-33</sup>

#### 4.2.1 Tube coagulase test direct from colonies

- use commercially available plasma and dilute according to manufacturer's instructions unless an alternative method has been validated
- label the test tubes with the organism to be tested as well as the control organisms
- emulsify representative colony/colonies of the test organism in the plasma. Incubate at 35-37°C and examine hourly up to 4hr. Do not shake or agitate the tube
- gently slant and examine for a clot which gels the whole contents of the tube or forms a loose web of fibrin

- if negative by the end of 4hr, incubate overnight at room temperature (22°C) and re-examine at 24hr. This is because a small proportion of strains require longer than 4hr for clot formation<sup>2</sup>

#### 4.2.2 Tube coagulase test direct from blood culture broth

This can be used for a rapid presumptive identification of *S. aureus* in blood culture broths.

- add 1 – 2 drops of the positive blood culture broth to 2mL of the diluted plasma in a tube or bijou
- incubate at 35-37°C and examine hourly up to 4hr
- examine for a clot which gels the whole contents of the tube or forms a loose web of fibrin
- if negative by the end of 4hr, incubate overnight at room temperature (22°C) and re-examine at 24hr. This is because a small proportion of strains require longer than 4hr for clot formation<sup>2</sup>

**Note:** Always check identification the following day from the culture plate using a slide coagulase or latex test.

##### Positive result

Formation of a clot up to 4hr at 37°C or following overnight incubation at room temperature (22°C).

##### Negative result


No clot, plasma moves freely at 4hr and 24hr incubation.

Table 1: Summary of coagulase test of some *Staphylococcus* species that may cause infections in humans<sup>1, 2</sup>.

Species	Tube coagulase test	Slide coagulase test
<i>Staphylococcus aureus</i> subspecies <i>aureus</i>	+	+
<i>Staphylococcus aureus</i> subspecies <i>anaerobius</i>	+	-
<i>Staphylococcus epidermidis</i>	-	-
<i>Staphylococcus haemolyticus</i>	-	-
<i>Staphylococcus saprophyticus</i> subspecies <i>saprophyticus</i>	-	-
<i>Staphylococcus schleiferi</i> subspecies <i>coagulans</i>	-	+
<i>Staphylococcus lugdunensis</i>	-	+
<i>Staphylococcus schleiferi</i> subspecies <i>schleiferi</i>	-/+	+
<i>Staphylococcus delphini</i> *	+	-

<i>Staphylococcus intermedius</i> *	+	v
<i>Staphylococcus hyicus</i> *	V	-
V= variable reaction                      - = negative reaction *rare clinical isolates                    +=positive reaction		

## Appendix: Coagulase test



### Note:

**Positive control:** *Staphylococcus aureus* NCTC 6571 or NCTC 12973

**Negative control:** *Staphylococcus haemolyticus* NCTC 11042

\* If blood cultures are tested directly by tube coagulase test and are found to be negative, follow the procedure for culture above for re-testing.

The flowchart is for guidance only.

## References

---

### Modified GRADE table used by UK SMI's when assessing references

Grading of Recommendations, Assessment, Development, and Evaluation (GRADE) is a systematic approach to assessing references. A modified GRADE method is used in UK SMI's for appraising references for inclusion. Each reference is assessed and allocated a grade for strength of recommendation (A-D) and quality of the underlying evidence (I-VI). A summary table which defines the grade is listed below and should be used in conjunction with the reference list.

Strength of recommendation		Quality of evidence	
A	Strongly recommended	I	Evidence from randomised controlled trials, meta-analysis and systematic reviews
B	Recommended but other alternatives may be acceptable	II	Evidence from non-randomised studies
C	Weakly recommended: seek alternatives	III	Non-analytical studies, for example, case reports, reviews, case series
D	Never recommended	IV	Expert opinion and wide acceptance as good practice but with no study evidence
		V	Required by legislation, code of practice or national standard
		VI	Letter or other

1. MacFaddin JF. Coagulase Test. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 105-19. **B, III**
2. Brown DF, Edwards Di, Hawkey PM, Morrison D, Ridgway GL, Towner KJ et al. Guidelines for the laboratory diagnosis and susceptibility testing of methicillin-resistant Staphylococcus aureus (MRSA). JAntimicrobChemother 2005;56:1000-18. **A, IV**
3. Williams REO HG. Determination of Coagulase and Alpha Haemolysin Production by Staphylococci. Br J Exp Pathol 1946;27:72-81. **B, I**
4. Katz DS. Coagulase Test Protocol. American Society for Microbiology peer-reviewed. 2012. **B, VIII**
5. Yrios JW. Comparison of rabbit and pig plasma in the tube coagulase test. JClinMicrobiol 1977;5:221-4. **B, II**
6. Sperber WH, Tatini SR. Interpretation of the tube coagulase test for identification of Staphylococcus aureus. ApplMicrobiol 1975;29:502-5. **B, II**

7. Cooke RP, Jenkins CT. Comparison of commercial slide agglutination kits with a tube coagulase test for the rapid identification of *Staphylococcus aureus* from blood culture. *JClinPathol* 1997;50:164-6. **B, II**
8. Ozen NS, Ogunc D, Mutlu D, Ongut G, Baysan BO, Gunseren F. Comparison of four methods for rapid identification of *Staphylococcus aureus* directly from BACTEC 9240 blood culture system. *Indian JMed Microbiol* 2011;29:42-6. **B, II**
9. Rossney AS, English LF, Keane CT. Coagulase testing compared with commercial kits for routinely identifying *Staphylococcus aureus*. *JClinPathol* 1990;43:246-52. **B, II**
10. Davis TE, Fuller DD, Aeschleman EC. Rapid, direct identification of *Staphylococcus aureus* and *Streptococcus pneumoniae* from blood cultures using commercial immunologic kits and modified conventional tests. *DiagnMicrobiolInfectDis* 1992;15:295-300. **B, II**
11. McDonald CL, Chapin K. Rapid identification of *Staphylococcus aureus* from blood culture bottles by a classic 2-hour tube coagulase test. *JClinMicrobiol* 1995;33:50-2. **B, II**
12. Kearns AM, Ganner M, Holmes A. The 'Oxford *Staphylococcus*': a note of caution. *J Antimicrob Chemother* 2006;58:480-1. **B, VII**
13. European Parliament. UK Standards for Microbiology Investigations (UK SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU *in vitro* Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes". 1998. **A, V**
14. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices 1998. 1-37. **A, V**
15. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 2009. **A, V**
16. Department for Transport. Transport of Infectious Substances, 2011 Revision 5. 2011. **A, V**
17. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2017-2018. 2017. **A, V**
18. Home Office. Anti-terrorism, Crime and Security Act. 2001. **A, V**
19. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive 2013. 1-35. **A, V**
20. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office 2003. **A, V**
21. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive 2005. **A, V**
22. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive 2008. **A, V**



23. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. MMWR Surveill Summ 2012;61:1-102. **B, IV**
24. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002 (as amended). HSE Books,. 2013. **A, V**
25. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books,. 2002. **A, V**
26. Health and Safety Executive. A Guide to Risk Assessment Requirements: Common Provisions in Health and Safety Law. HSE Books,. 2002. **A, V**
27. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books 2003. **A, V**
28. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets 2000. **A, V**
29. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 2005. 1-14. **A, V**
30. Department of Health. Transport of Infectious Substances. Best Practice Guidance for Microbiology Laboratories. Department of Health. 1-13. 2007. **A, V**
31. Barrow GI, Feltham RKA. Bacterial characters and characterization. Cowan and Steel's Manual for the Identification of Medical Bacteria. 3rd ed. Cambridge: Cambridge University Press; 2003. p. 31-2, 226. **B, III**
32. Evans JB, Grun, L. . Subcommittee on taxonomy of staphylococci and micrococci - minutes of first meeting. International bulletin of bacteriological nomenclature and taxonomy 1965;15:108. **C, III**
33. Zarzour JY, Belle EA. Evaluation of three test procedures for identification of Staphylococcus aureus from clinical sources. JClinMicrobiol 1978;7:133-6. **C, III**