



UK Health
Security
Agency

UK Standards for Microbiology Investigations

Identification of *Neisseria* species



National Institute for Health and Care Excellence (NICE) has renewed accreditation of the process used by the UK Health Security Agency to produce UK Standards for Microbiology Investigations (UK SMIs). The renewed accreditation is valid until 30 June 2026 and applies to guidance produced using the processes described in 'UK Standards for Microbiology Investigations Development Process' (2021). The original accreditation term began on 1 July 2011.

Consultation between 7 September 2023 to 22 September 2023

Acknowledgments

UK Standards for Microbiology Investigations (UK SMIs) are developed under the auspices of UKHSA working in partnership with the partner organisations whose logos are displayed below and listed on [the UK SMI website](#). UK SMIs are developed, reviewed and revised by various working groups which are overseen by a steering committee (see [the Steering Committee page on GOV.UK](#)).

The contributions of many individuals in clinical, specialist and reference laboratories who have provided information and comments during the development of this document are acknowledged. We are grateful to the medical editors for editing the medical content.

UK SMIs are produced in association with:



Displayed logos correct as of March 2021

UKHSA publications gateway number:

Contents

Acknowledgments	2
Contents	3
Amendment table	4
1 General information	5
2 Scientific information	5
3 Scope of document	5
4 Introduction	5
5 Technical information and limitations	6
6 Safety considerations	6
7 Target organisms	7
8 Identification	7
9 Reporting	13
10 Referral to reference laboratories	15
Algorithm 1: Identification of <i>Neisseria gonorrhoea</i>	16
Algorithm 2: Identification of non-gonococcal <i>Neisseria</i> species	17
References	18

Consultation between 7 September 2023 to 22 September 2023

Amendment table

Each UK SMI document has an individual record of amendments. The amendments are listed on this page. The amendment history is available from standards@ukhsa.gov.uk.

Any alterations to this document should be controlled in accordance with the local document control process.

Amendment number/date	x/dd.mm.yy
Issue number discarded	
Insert issue number	
Anticipated next review date*	dd.mm.yy
Section(s) involved	Amendment

Consultation between 7 September 2023 to 22 September 2023

1 General information

[View general information](#) related to UK SMIs.

2 Scientific information

[View scientific information](#) related to UK SMIs.

3 Scope of document

This UK Standards for Microbiology Investigations (UK SMI) document describes the identification of *Neisseria* species and includes routine culture, microscopy, oxidase test and MALDI-TOF MS for identification. It also covers biochemical tests and molecular methods for confirmation.

This document describes the differentiation of pathogenic *Neisseria* species from non-pathogenic *Neisseria* species and the related genera of *Moraxella* and *Kingella*. The identification of these genera is covered in [ID 11 - Identification of *Moraxella* species and morphologically similar organisms](#) and [ID 12 – Identification of *Haemophilus* species and the HACEK group of organisms](#).

This document does not focus on the screening of *Neisseria gonorrhoea* and *Neisseria meningitidis* or antimicrobial susceptibility testing of *Neisseria* species. The screening of *N. meningitidis* is covered in [B 51 - Screening for *Neisseria meningitidis*](#).

Some of the *Neisseria* species have been reclassified, and the updated nomenclature of these species have been included in this document for reference.

This UK SMI should be used in conjunction with other UK SMIs.

4 Introduction

4.1 Taxonomy and characteristics

The genus *Neisseria* comprises gram-negative bacteria belonging to the family *Neisseriaceae*, order *Neisseriales* within the phylum β -Proteobacteria (1-4). There are currently more than 30 *Neisseria* species and 3 subspecies of which may be isolated from humans and animals (1,2). The following species: *Neisseria ovis*, *Neisseria cuniculi* and *Neisseria caviae* have been reclassified to *Moraxella ovis*, *Moraxella cuniculi* and *Moraxella caviae*, respectively (1).

The clinically important species are *Neisseria gonorrhoeae* (*N. gonorrhoeae*), and *Neisseria meningitidis* (*N. meningitidis*). The respective species are closely related but cause entirely different diseases with distinct clinical pathologies. *N. gonorrhoeae* is an obligate pathogen that causes the sexually transmitted infection gonorrhoea.

N. meningitidis is an opportunistic pathogen that colonises the nasopharyngeal mucosa and has the potential to cause meningococcal disease which includes meningitis and septicemia.

The other *Neisseria* species such as *Neisseria lactamica* and *Neisseria cinerea* are generally considered commensals, but have been implicated as causes of infection in patients who are immunocompromised (5). More recent *Neisseria* species implicated in human disease include *N. brasiliensis*, *N. dumasiana*, *N. oralis*, *N. shayegani*, *N. wadsworthii* and *N. skkuensis* (1).

Characteristics

Neisseria species are Gram negative cocci, 0.6 to 1.0µm in diameter, occurring singly but more often in pairs with adjacent sides flattened; except *Neisseria elongata*, *Neisseria weaver*, *Neisseria bacilliformis* and *Neisseria shayegani* that consist of rods, 0.5µm wide, often arranged as diplococci or in short chains (6-9). They are non-motile (10). Some species produce a greenish-yellow carotenoid pigment, and some may be nutritionally fastidious and haemolytic. Some species are saccharolytic. The optimum growth temperature is 35 to 37°C. *Neisseria* are oxidase positive and catalase positive (except *Neisseria elongata*).

5 Technical information and limitations

The advancement in molecular typing revealed that *Neisseria* species are larger and more diverse than previously thought (1,10). This led to the discovery of many novel species and the reclassification and nomenclature changes of others (1).

The changes made in the taxonomy of the *Neisseria* genus need to be reflected in the databases of the identification tools used in laboratories, this is particularly important for species that are closely related to *N. gonorrhoea* and *N. meningitidis* as a lack of match of these closely related species to the database and subsequent identification can lead to uncertainty and misidentification resulting in serious consequences (11-13).

Note: The social consequences to the patient and the organization of an incorrect diagnosis of gonorrhoea disease as a result of misidentification should not be underestimated.

6 Safety considerations

The section covers specific safety considerations (14-35) related to this UK SMI, and should be read in conjunction with the general [safety considerations on GOV.UK](#).

N. meningitidis is a Hazard group 2 organism, the processing of diagnostic samples should be carried out at Containment Level 2.

Due to the severity of the disease and the risks associated with generating aerosols, any manipulation of suspected isolates of *N. meningitidis* should always be undertaken in a microbiological safety cabinet until *N. meningitidis* has been ruled out (as must any laboratory procedure giving rise to infectious aerosols) (31).

N. meningitidis causes severe and sometimes fatal disease. Laboratory acquired infections have been reported (36,37). The organism infects primarily by the respiratory route. An effective vaccine is available for most meningococcal groups. Vaccination is required for laboratory staff routinely working with the organism.

N. gonorrhoeae is also a Hazard group 2 organism which is responsible for the sexually transmitted infection, gonorrhoea and can also cause eye or throat infection - which is the most likely risk to laboratory workers through either vertical transmission, poor laboratory practice or inhalation of aerosols.

Refer to current guidance on the safe handling of all organisms documented in this UK SMI.

The above guidance should be supplemented with local COSHH and risk assessments.

Compliance with postal and transport regulations is essential.

7 Target organisms

Please refer to Table 1 for *Neisseria* species that have been associated with human disease (1).

Other organisms which may be misidentified as *Neisseria* species are *Moraxella catarrhalis* and *Kingella denitrificans* (38,39).

8 Identification

Culture-based methods remain the gold standard for identification, with the integration of faster identification techniques such as MALDI-TOF MS improving the accuracy of identification. There is also a growing shift towards molecular methods for identification. However, these techniques require specialised laboratories, trained staff and expensive reagents which may not be available to all routine laboratories.

8.1 Culture methods

Culture methods provide presumptive identification of *Neisseria* species based on colony morphology (in some cases - Gram stain) and oxidase followed by identification via MALDI-TOF MS. Additional confirmation can be performed using biochemical or molecular tests.

In cases where confirmation is not possible and further identification is required, isolates should be referred to the appropriate reference laboratory.

8.1.1 Bacterial growth media

Some *Neisseria* species including *N. gonorrhoeae* and *N. meningitidis* are fastidious and require enriched media for growth. They grow best in aerobic conditions at temperatures of 35 to 37°C with 5 to 10% CO₂ (40). Colonies usually appear within 18 to 48 hours of incubation and vary in morphological appearance.

8.1.1.1 Primary agar

Whole blood agar or heated blood (chocolate) incubated for 18 to 48 hours in 5 to 10% CO₂ at 35 to 37°C (12). The media usually consist of Columbia agar base supplemented with 5% horse blood or chocolatised horse blood.

Note: *N. gonorrhoeae* grows poorly on blood agar, so a non-selective GC agar with lysed or chocolatised horse blood should be used instead (41).

8.1.1.2 Selective agar

GC selective agar incubated for >40 hours in 5 to 10% CO₂ at 35 to 37°C. This selective agar is primarily used for the selective isolation of *N. gonorrhoeae* but can also be used for the isolation of *N. meningitidis*.

Note: *Neisseria* species except *Neisseria lactamica* generally do not grow well on Thayer-Martin based GC selective agar and can be differentiated from *N. gonorrhoeae* and *N. meningitidis* using methods such as MALDI-TOF MS or biochemical tests.

8.1.2 Colonial appearance

Neisseria species are usually pigmented and opaque. However, both *N. gonorrhoeae* and *N. meningitidis* form smooth, round, moist, uniform grey/brown colonies with a greenish colour underneath on primary isolation medium. Table 3 details the colony morphology of *Neisseria* species.

8.2 Microscopic appearance

8.1.2.1 Gram stain

Please refer to 2.39 – Staining procedures

Neisseria species

Gram-negative cocci arranged in pairs with long axes parallel or gram-negative rods that are arranged in chains or as diplococci, table 3 details the microscopic appearance of *Neisseria* species in Gram's stain.

Note: Gram-stain is often omitted from the identification process if isolates are going to be identified using MALDI-TOF MS.

Table 1. Microscopic and colonial morphology of *Neisseria* species (1)

Please note that the information in this table provides general characteristics of colony appearance and can vary among different strains and culture conditions.

Species	Colonies	Additional comments
<i>Neisseria gonorrhoeae</i>	Diplococci with concave adjacent sides. Smooth, round, moist, uniform grey/brown with a greenish colour underneath	Non-haemolytic. No pigmentation. Poor growth on blood agar when the medium is very fresh, or the number of bacteria present in the sample is especially high. Autolysis and sticky colonies with prolonged growth.
<i>Neisseria meningitidis</i>	Diplococci. Similar to <i>N. gonorrhoeae</i>	Non-haemolytic on blood agar. No pigmentation. Autolysis with prolonged growth
<i>Neisseria lactamica</i>	Diplococci. Colonies less moist and smaller than <i>N. gonorrhoeae</i> and <i>N. meningitidis</i>	Haemolytic on horse blood agar*. Yellow pigmentation
<i>Neisseria cinerea</i>	Diplococci/scattered clusters. Small, greyish white with entire edges, and slightly granular	Non-haemolytic. Yellow pigmentation*
<i>Neisseria elongata</i>	Small slender rods that occur in chains. Small, greyish white, shiny opaque colonies, low-hemispherical with an entire edge	Non-haemolytic with some pitting of the agar. Yellow pigmentation*
<i>Neisseria elongata</i> sub. <i>elongata</i>	Flat colonies	Non-haemolytic
<i>Neisseria elongata</i> subsp. <i>glycolytica</i>	Similar to <i>N. elongata</i> colonies.	Haemolysis varies. Yellow pigmentation*. Relatively large grey, opaque, moderately raised with flat top and smooth with a soft homogenous consistency on blood agar.
<i>Neisseria elongata</i> subsp. <i>nitroreducens</i>	Similar to <i>N. elongata</i> colonies.	
<i>Neisseria sicca</i>	Cocci occurring in pairs and tetrads. small round colonies, having a smooth surface and an entire edge	Haemolytic*. Yellow pigmentation*. Colonies increase in size, and appear raised, rough, and black after 24 hrs. Very firm to the medium.
<i>Neisseria mucosa</i>	Diplococci. Large, mucoid, and often adherent	No haemolysis*. No pigmentation or greyish to buff yellow.

Consultation between 7 September 2023 to 22 September 2023

<i>Neisseria canis</i>	Diplococci/rarely in tetrads. Smooth, butyrous with a light-yellow tinge	No haemolysis. No pigmentation.
<i>Neisseria flava</i>	Diplococci. Discrete, opaque, pale-yellow, slightly flatter than <i>N. meningitidis</i> colonies.	Yellow pigmentation*
<i>Neisseria subflava</i>	Cocci occurring in pairs and tetrads. Smooth, transparent, or opaque, often adherent	No haemolysis. Yellow pigmentation. They tend to resist Gram decolourisation
<i>Neisseria bacilliformis</i>	Small rods. Round, smooth, glistening, light grey	
<i>Neisseria weaveri</i>	Broad, plump, medium-to-large, straight rods of varying length in chains or longer rods. Smooth, flat, somewhat glistening with an entire border.	Haemolytic. Colonies are variable in size and increase after 24hrs
<i>Neisseria flavescens</i>	Cocci occurring in pairs or tetrads. Smooth and opaque.	Non-haemolytic. Golden/yellow pigmentation.
<i>Neisseria oralis</i>	Cocci occurring in chains. Small, circular, entire, raised, moist, and yellow.	Weakly haemolytic.
<i>Neisseria shayegani</i>	Rod-shaped and long. Small, circular, entire, convex, moist, light yellow/grey.	Non-haemolytic
<i>Neisseria wadsworthii</i>	Diplococci in chains. Small, circular, entire, convex, moist, light yellow/orange	Non-haemolytic
<i>Neisseria zoodegmatidis</i>	Coccoid rods. Circular, convex, entire, opaque, shiny, and smooth	Haemolytic. No pigmentation
<i>Neisseria animalis</i>	Coccoid rods. Colonies are circular, convex, entire, opaque, shiny, and smooth	Haemolytic. No pigmentation.
<i>Neisseria dumasiana</i>	Coccoid to coccobacilli, may be present in pairs. Grey, moist, circular, convex, entire	Non-haemolytic. Grey pigmentation
<i>Neisseria brasiliensis</i>	Diplococci, brownish colonies	
<i>Neisseria Skkuensis</i>	Small, round, and light grey	
<i>Neisseria polysaccharea</i>	Cocci arranged in pairs or tetrads. Relatively small (2mm) yellowish colonies	Non-haemolytic. Large amounts of polysaccharides produced

Consultation between 7 September 2022 to 22 September 2023

<i>Neisseria caviae</i> *	Diplococci with adjacent sides flattened. Small (2mm), circular, convex with entire edge, and a smooth glistening surface. butyrous becoming viscid.	Weakly haemolytic*. Light caramel-light brown pigmentation
<i>Neisseria cuniculi</i> *	Oval cocci, small, and smooth	Haemolytic
<i>Neisseria ovis</i> *	Diplococci. Grey, opaque, convex	Haemolytic
* Reclassified		

8.3 Oxidase test

Please refer to [TP 26 – Oxidase test](#)

Neisseria species are oxidase positive.

Note: *Kingella* species and *M. catarrhalis* are also oxidase positive and can be misidentified as *Neisseria* (38,39).

8.4 Matrix-assisted laser desorption/ionisation - time of flight mass spectrometry (MALDI-TOF MS)

MALDI-TOF MS is currently used as the primary method for the identification of *Neisseria* species, while biochemical tests and molecular methods serve as an alternative identification approach for centres without MALDI-TOF MS or are employed as a confirmatory method (41).

Although the problem of the *Neisseria* genus study is complex, MALDI-TOF MS has been developed and validated to differentiate the clinically important species, *N. gonorrhoeae* and *N. meningitidis*.

MALDI-TOF MS has excellent performance for *N. gonorrhoeae* identification (12,42) and is highly accurate for the identification of *N. meningitidis* however closely related *Neisseria* species such as *Neisseria polysaccharea* and *Neisseria cinerea* may be misidentified as *N. meningitidis* (11,12,43). While the identification of non-pathogenic *Neisseria* to species level is generally not required, the misidentification of commensal strains as *N. meningitidis* can have serious health and social consequences.

Therefore, confirmatory testing with biochemical tests or molecular methods may be required if a consistent identification is not achieved by MALDI-TOF MS (43).

Continual improvement of MALDI-TOF MS requires enriching its database with spectra from closely related and poorly represented *Neisseria* species (11). This is essential to increase specificity for *N. meningitidis* (11). The ongoing efforts to expand

and refine the database will ultimately improve reliability and accuracy of MALDI-TOF MS for identification.

8.5 Further identification

8.5.1 Biochemical tests and commercial identification systems

Biochemical tests including commercial identification kits provide basic biochemical information that can aid in the identification of *Neisseria* species. However, relying solely on these tests is insufficient for accurate identification of *Neisseria* species. Therefore, these tests are not considered reliable for the primary identification of *Neisseria* species.

Refer to manufacturer's guidance or the Manual of Clinical Microbiology for the biochemical properties of *Neisseria* species (44).

Commercially available kits can be used for confirmation of MALDI-TOF MS results. The accuracy of these kits has not been fully determined for species other than *N. gonorrhoeae* and *N. meningitidis* therefore, all results obtained should be interpreted with caution.

Laboratories should follow manufacturers' instructions and rapid tests and kits should be validated and be shown to be fit for purpose prior to use.

Currently, there are limited immunological kits available for the identification of *Neisseria* species. However, biochemical kits can be used as an alternative. Many of the biochemical kits allow for the combined detection of preformed enzymes and carbohydrate utilisation.

In addition, commercial latex or slide agglutination kits can be used for further characterisation of *N. meningitidis* to serogroup level (45,46). The latex agglutination kits are designed for direct use on CSF or serum but will also work for cultures. Slide agglutinating sera are for use on cultures only. Heated clinical samples or formalin treated suspensions of cultures should be processed within microbiological safety cabinets to reduce aerosols.

8.5.2 Molecular methods

Molecular techniques have made identification of many species more rapid and precise than is possible with phenotypic techniques. The routine implementation of molecular methods can be challenging, as not all clinical laboratories have access to the different molecular methods. Therefore, in such cases significant isolates identified by MALDI-TOF MS should be sent to appropriate reference laboratories for further testing and confirmation of results if required.

8.3.2.1 Polymerase Chain Reaction (PCR)

PCR is mainly used as complementary or confirmatory testing method for the identification of *Neisseria* species following MALDI-TOF MS results (47). Laboratories that are unable to perform PCR for confirmation of MALDI-TOF MS results can carry

out immunological or biochemical testing using available commercial kits and/or send isolates to reference laboratories for further testing as required.

8.3.2.2 Next generation sequencing (NGS)

With the increased availability of NGS technologies, there may be a shift towards their utilisation for the identification of *Neisseria* species alongside other target pathogens in the future. However currently these technologies are largely restricted to reference units.

Whole genome sequencing (WGS) is routinely used by UKHSA and has greatly improved surveillance capabilities and monitoring trends in antimicrobial resistance. WGS has replaced traditional phenotypic and polymerase chain reaction (PCR) methods for routine surveillance. It also has high discriminatory power and can provide in-depth genetic analysis and identification (10). Therefore, it has the potential to be an alternative to techniques like MALDI-TOF MS for the identification of *Neisseria* including unknown *Neisseria* species (10).

8.4 Storage and referral

Short term storage – isolates should be kept in a viable state on heated blood (chocolate) agar slopes.

Long term storage – isolates should be frozen at -20°C to -80°C in glycerol based medium or cry-beads (41).

Note: *N. gonorrhoeae* storage is recommended at -70°C

9 Reporting

9.1 Infection Specialist

Inform the infection specialist of all confirmed *N. meningitidis* isolates, and of all *Neisseria* species isolated from normally sterile sites, or in cases of invasive infection.

The infection specialist should also be informed if the request bears relevant information, for example:

- cases of meningitis, septicaemia (especially with purpuric rash)
- investigation of *N. meningitidis* outbreak, or of the carrier state

Inform the infection specialist of all confirmed *N. gonorrhoeae* isolates, and of all *Neisseria* species from:

- minors
- cases of sexual assault, rape, or abuse

- cases of *N. gonorrhoeae* isolated from normally sterile sites or from invasive infection – also send to the appropriate reference laboratory
- Multi-drug resistant isolates of GC from all sites.

Follow local protocols for reporting to clinician.

9.2 Presumptive identification

It is not recommended that presumptive identifications for *Neisseria meningitidis* or *Neisseria gonorrhoea* are reported due to the risk associated with mis-diagnosing both infections before full identification are obtained by MALDI-TOF MS or PCR. For centres where these tests are not available a presumptive identification can be made using 2 or 3 biochemical and immunological tests or commercial kits but should be confirmed by accurate methods.

9.3 Confirmation of identification

Any one of the approaches listed below can be taken to confirm the identity of the *Neisseria* species following identification processes as outlined in this document and/or Reference Laboratory report

1. MALDI-TOF MS confirmation
2. Molecular confirmation using PCR

Note: Commercial latex kit or slide agglutination reagent is an additional step for further confirmation and characterisation of *N. meningitidis* to serogroup

For confirmation and identification refer to section 10.

9.4 Health Protection Team (HPT)

Refer to local agreements in devolved administrations.

9.5 UK Health Security Agency

Refer to current guidelines on Second Generation Surveillance System (SGSS) reporting (30).

9.6 Infection prevention and control team

Inform the infection prevention and control team of presumptive and confirmed isolates of *N. meningitidis*.

10 Referral to reference laboratories

For information on the tests offered, turnaround times, transport procedure and the other requirements of the reference laboratory [see user manuals and request forms](#)

Contact appropriate reference laboratory for information on the tests available, turnaround times, transport procedure and any other requirements for sample submission:

[England](#)

[Wales](#)

[Scotland](#)


[Northern Ireland](#)

Note: In case of sending away to laboratories for processing, ensure that specimen is placed in appropriate package and transported accordingly.

Consultation between 7 September 2023 to 22 September 2023

Algorithm 1: Identification of *Neisseria gonorrhoea*


An accessible text description of this flowchart is provided with this document



The flowchart is for guidance only.

Algorithm 2: Identification of non-gonococcal *Neisseria* species

An accessible text description of this flowchart is provided with this document



The flowchart is for guidance only.

References

An explanation of the reference assessment used is available in the [scientific information section on the UK SMI website](#).

1. Parte AC and others. 'List of Prokaryotic names with Standing in Nomenclature (LPSN) moves to the DSMZ'. *International Journal of Systematic and Evolutionary Microbiology* 2020: volume 70, issue 11, pages 5607-12. <https://doi.org/10.1099/ijsem.0.004332>
2. Schoch CL and others. 'NCBI Taxonomy: a comprehensive update on curation resources and tools'. *Database (Oxford)* 2020: volume 2020. 10.1093/database/baaa062
3. Hung MC, Christodoulides M 'The biology of *Neisseria adhesins*'. *Biology (Basel)* 2013: volume 2, issue 3, pages 1054-109. 10.3390/biology2031054
4. Jolley KA and others. 'Open-access bacterial population genomics: BIGSdb software, the PubMLST.org website and their applications'. *Wellcome Open Research* 2018: volume 3, pages 124. 10.12688/wellcomeopenres.14826.1
5. Dorey RB and others. 'The nonpathogenic commensal *Neisseria*: friends and foes in infectious disease'. *Current Opinion in Infectious Diseases* 2019: volume 32, issue 5, pages 490-6. 10.1097/qco.0000000000000585
6. Bergey DH and Holt JG 'Bergey's Manual of Determinative Bacteriology'. 9th edition. Baltimore: Williams and Wilkins, 1994. pages. 34-540
7. Han XY and others. '*Neisseria bronchialis* sp. nov. isolated from human infections'. *Journal of Clinical Microbiology* 2006: volume 44, issue 2, pages 474-9. 10.1128/JCM.44.2.474-479.2006
8. Andersen BM and others. '*Neisseria weaveri* sp. nov., formerly CDC group M-5, a gram-negative bacterium associated with dog bite wounds'. *Journal of Clinical Microbiology* 1993: volume 31, issue 9, pages 2456-66
9. Wolfgang W and others. '*Neisseria wadsworthii* sp. nov. and *Neisseria shayegani* sp. nov., isolated from clinical specimens'. *International Journal of Systematic and Evolutionary Microbiology* 2011: volume 61, issue Pt 1, pages 91-8. 10.1099/ij.s.0.022426-0
10. Danler CM 'Neisseria species and their complicated relationships with human health'. *Microbiology Australia* 2021: volume 42, issue 2, pages 79-83. <https://doi.org/10.1071/MA21024>
11. Hong E and others. 'Identification of *Neisseria meningitidis* by MALDI-TOF MS may not be reliable'. *Clinical Microbiology and Infection* 2019: volume 25, issue 6, pages 717-22. 10.1016/j.cmi.2018.09.015

12. Morel F and others. 'Use of Andromas and Bruker MALDI-TOF MS in the identification of *Neisseria*'. *European Journal of Clinical Microbiology and Infectious Disease* 2018: volume 37, issue 12, pages 2273-7. 10.1007/s10096-018-3368-6
13. Schweitzer VA and others. 'Identification of *Neisseria gonorrhoeae* by the Bruker Biotyper Matrix-Assisted Laser Desorption Ionization-Time of Flight Mass Spectrometry System Is Improved by a Database Extension'. *Journal of Clinical Microbiology* 2016: volume 54, issue 4, pages 1130-2. 10.1128/jcm.00016-16
14. Advisory Committee on Dangerous Pathogens. '[The Approved List of Biological Agents](#)'. Health and Safety Executive 2021. pages 1-39 ++
15. British Standards Institution (BSI). '[BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets](#)'. 2000 ++
16. British Standards Institution (BSI). '[BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets](#)'. Recommendations and guidance. 2005. pages 1-14 ++
17. Centers for Disease Control and Prevention. '[Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories](#)'. *The Morbidity and Mortality Weekly Report (MMWR) Surveillance Summaries* 2012, volume 61, pages 1-102 +
18. Department for Transport and others. '[Transport of infectious substances UN2814, UN2900 and UN3373 Guidance note number 12012 \(revision 7\)](#)'. 2013 ++
19. Department of Health. '[Health Protection Legislation \(England\) Guidance](#)'. pages 1-112. 2010 ++
20. Gizzie N, Adukwu E 'Evaluation of Liquid-Based Swab Transport Systems against the New Approved CLSI M40 A2 Standard'. *Journal of Clinical Microbiology* 2016: volume 54, issue 4, pages 1152-6. 10.1128/JCM.03337-15 +
21. Health and Safety Executive. '[Managing risks and risk assessment at work](#)' (accessed 28/07/2021) ++
22. Health and Safety Executive. '[Safe use of pneumatic air tube transport systems for pathology specimens](#)'. 2009 ++
23. Health and Safety Executive. '[Control of Substances Hazardous to Health Regulations](#)'. *The Control of Substances Hazardous to Health Regulations 2002 (as amended)*. Approved Code of Practice and guidance L5 (sixth edition). Health and Safety Executive (HSE) Books. 2013 ++
24. Health and Safety Executive. '[Risk assessment: A brief guide to controlling risks in the workplace](#)'. Health and Safety Executive (HSE). 2014 ++
25. Health and Safety Executive. Advisory Committee on Dangerous Pathogens. '[Management and operation of microbiological containment laboratories](#)'. Health and Safety Executive (HSE). 2019 ++

26. Health Services Advisory Committee. '[Safe working and the prevention of infection in clinical laboratories and similar facilities](#)'. Health and Safety Executive (HSE) Books. 2003 ++
27. Home Office. '[Public Health Act \(Northern Ireland\)](#)'. 1967 Chapter 36 ++
28. Home Office. '[Anti-terrorism, Crime and Security Act](#)'. 2001 ++
29. Official Journal of the European Communities. 'Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro diagnostic medical devices'. 1998. pages 1-37 ++
30. Public Health England. '[Laboratory reporting to Public Health England: a guide for diagnostic laboratories](#)'. Public Health of England (PHE). 2020. pages 1-31 ++
31. Scottish Government. '[Public Health \(Scotland\) Act](#)'. 2008 ++
32. The Royal College of Pathologists. '[The retention and storage of pathological records and specimens \(5th edition\)](#)'. pages 1-59. 2015 ++
33. The Welsh Assembly Government. '[Health Protection Legislation \(Wales\) Guidance](#)'. 2010 ++
34. Tyrrell KL and others. 'Comparison of the Copan eswab System with an Agar Swab Transport System for Maintenance of Fastidious Anaerobic Bacterium Viability'. *Journal of Clinical Microbiology* 2016: volume 54, issue 5, pages 1364-7. 2+
10.1128/JCM.03246-15
35. World Health Organization. '[Guidance on regulations for the transport of infectious substances 2019-2020](#)'. World Health Organisation (WHO). 2019 ++
36. Sejvar JJ and others. 'Assessing the risk of laboratory-acquired meningococcal disease'. *Journal of Clinical Microbiology* 2005: volume 43, issue 9, pages 4811-4. 10.1128/JCM.43.9.4811-4814.2005
37. Bhatti AR and others. 'A laboratory-acquired infection with *Neisseria meningitidis*'. *Journal of Infection* 2002: volume 4, issue 3, pages 247-52
38. Henriksen SO, Bøvre K 'The Taxonomy of the Genera *Moraxella* and *Neisseria*'. *Microbiology* 1968: volume 51, issue 3, pages 387-92.
<https://doi.org/10.1099/00221287-51-3-387>
39. Smith JJS, Lapage SP 'Transfer of Some Saccharolytic *Moraxella* Species to *Kingella* Henriksen and Bøvre 1976, with Descriptions of *Kingella indologenes* sp. nov and *Kingella denitrificans* sp. nov'. *International Journal of Systematic and Evolutionary Microbiology* 1976: volume 26, issue 4, pages 451-8.
<https://doi.org/10.1099/00207713-26-4-451>
40. Mehrabi Tavana A and others. 'Characteristics of *Neisseria* Species Colonized in the Human's Nasopharynx'. *Jundishapur Journal of Microbiology* 2020: volume 13, issue 6, pages e99915. 10.5812/jjm.99915

41. Villa L and others. 'Long term storage of fastidious bacteria (*Neisseria* spp. and *Haemophilus* spp.) with swab preservation at -80°C '. *Journal of Microbiological Methods* 2020: volume 175, pages 105969.

<https://doi.org/10.1016/j.mimet.2020.105969>

42. Buchanan R and others. 'Matrix-assisted laser desorption-ionization time-of-flight mass spectrometry for the identification of *Neisseria gonorrhoeae*'. *Clinical Microbiology and Infection* 2016: volume 22, issue 9, pages 815.e5-.e7.

<https://doi.org/10.1016/j.cmi.2016.06.010>

43. Unalan-Altintop T and others. 'A diagnostic challenge in clinical laboratory: Misidentification of *Neisseria subflava* as *Neisseria meningitidis* by MALDI-TOF MS'. *Acta Microbiologica et Immunologica Hungarica* 2020, pages 1-3.

10.1556/030.2020.01039

44. Karen CC and Michael AP 'Manual of Clinical Microbiology'. 2023

45. UK Health Security Agency (UKHSA). '[Meningococcal reference unit: user manual](#)'. 2023

46. Ali RH and others. 'Production and Characterization of Bacterial Ghost Vaccine against *Neisseria meningitidis*'. *Vaccines (Basel)* 2022: volume 11, issue 1. 10.3390/vaccines11010037

47. Public Health England (PHE). '[Guidance for the detection of gonorrhoea in England](#)'. 2021

Consultation between 7 September 2023 to 22 September 2023