




# UK Standards for Microbiology Investigations

## ONPG ( $\beta$ -Galactosidase) test



## Acknowledgments

UK Standards for Microbiology Investigations (UK SMIs) are developed under the auspices of UKHSA working in partnership with the partner organisations whose logos are displayed below and listed on [the UK SMI website](#). UK SMIs are developed, reviewed and revised by various working groups which are overseen by a [steering committee](#).

The contributions of many individuals in clinical, specialist and reference laboratories who have provided information and comments during the development of this document are acknowledged. We are grateful to the medical editors for editing the medical content.

UK SMIs are produced in association with:



Displayed logos correct as of December 2024

## Contents

<b>Acknowledgments</b> .....	<b>2</b>
<b>Contents</b> .....	<b>3</b>
<b>Amendment table</b> .....	<b>4</b>
<b>1 General information</b> .....	<b>6</b>
<b>2 Scientific information</b> .....	<b>6</b>
<b>3 Scope of document</b> .....	<b>6</b>
<b>4 Introduction</b> .....	<b>6</b>
<b>5 Technical information/limitations</b> .....	<b>7</b>
<b>6 Safety considerations</b> .....	<b>7</b>
<b>7 Reagents and equipment</b> .....	<b>7</b>
<b>8 Quality control organisms</b> .....	<b>8</b>
<b>9 Procedure and results</b> .....	<b>8</b>
<b>Algorithm: ONPG (<math>\beta</math>-Galactosidase) test</b> .....	<b>9</b>
<b>References</b> .....	<b>10</b>

## Amendment table

Each UK SMI document has an individual record of amendments. The amendments are listed on this page. The amendment history is available from [standards@ukhsa.gov.uk](mailto:standards@ukhsa.gov.uk).

Any alterations to this document should be controlled in accordance with the local document control process.

Amendment number/date	8/28.02.25
Issue number discarded	4
Insert issue number	4.1
<b>Section(s) involved</b>	<b>Amendment</b>
Whole document.	<p><b>This is an administrative point change.</b></p> <p><b>The content of this UK SMI document has not changed.</b></p> <p><b>The last scientific and clinical review was conducted on 03/12/2018</b></p> <p>Hyperlinks throughout document updated to Royal College of Pathologists website.</p> <p>Public Health England replaced with UK Health Security Agency throughout the document, including the updated Royal Coat of Arms</p> <p>Partner organisation logos updated.</p> <p>Broken links to devolved administrations replaced.</p> <p>References to NICE accreditation removed.</p> <p>Scope and Purpose replaced with General and Scientific information to align with current UK SMI template.</p>

Amendment number/date	7/03.12.18
Issue number discarded	3
Insert issue number	4
Anticipated next review date*	03.12.21
<b>Section(s) involved</b>	<b>Amendment</b>

Whole document.	Document updated. Technical limitations updated with subheadings. References updated with grades.
-----------------	---

\*Reviews can be extended up to 5 years where appropriate

# 1 General information

[View general information](#) related to UK SMIs.

# 2 Scientific information

[View scientific information](#) related to UK SMIs.


# 3 Scope of document

This document covers the procedure for ONPG test. The test is important in differentiating among the Enterobacteriaceae which are commonly classified according to their ability to ferment lactose<sup>1</sup>. It is also used to differentiate *Neisseria lactamica* from other fastidious *Neisseria* species.

This UK SMI should be used in conjunction with other UK SMIs.

# 4 Introduction

The ONPG (o-nitrophenyl- $\beta$ -D-galactopyranoside) test is used to determine the presence or absence of the enzyme  $\beta$ -galactosidase in an organism<sup>2</sup>. The presence of two enzymes, permease and  $\beta$ -galactosidase, are required to demonstrate lactose fermentation. Permease allows the lactose to enter the bacterial cell. In lactose-fermenting bacteria the breakdown of lactose to glucose and galactose involves the enzyme beta-galactosidase<sup>3</sup>. True lactose non-fermenters do not possess either of these enzymes. Late lactose fermenting organisms do not have permease, but do possess  $\beta$ -galactosidase. ONPG is similar in structure to lactose. If  $\beta$ -galactosidase is present, the colourless ONPG is split into galactose and o-nitrophenol, a yellow compound<sup>4,5</sup>. The reaction is shown as follows:



**Note:** “ONPG” (also known as “2-Nitrophenyl  $\beta$ -D-galactopyranoside”) is a Chemical analog of the sugar lactose and is hydrolysed by the enzyme lactase. Like  $\beta$ -galactosidase, lactase breaks lactose down into galactose and glucose.

## 5 Technical information/limitations

### 5.1 Growth media

The test should be performed, where possible, from a non-selective medium. If the test is performed from selective agar, a purity plate must be included to check for purity of the organism. Organisms that have grown on glucose containing media show less reactivity than those grown on lactose containing media. Glucose inhibits  $\beta$ -galactosidase.

### 5.2 Pigmentation in organisms

The test cannot be performed on organisms containing a yellow pigment or other coloured pigmentation as it makes it difficult to read the test<sup>5</sup>.

### 5.3 Interpretation of results

The ONPG solution must be correctly buffered to prevent false negative and false positive reactions.

A heavy inoculum is necessary to obtain a high concentration of enzyme.

Discard the substrate if it looks yellow prior to inoculation.

## 6 Safety considerations<sup>6-23</sup>

Refer to current guidance on the safe handling of all organisms and reagents documented in this UK SMI.

All work likely to generate aerosols must be performed in a microbiological safety cabinet.

The above guidance should be supplemented with local COSHH and risk assessments.

Compliance with postal and transport regulations is essential.

## 7 Reagents and equipment

Discrete bacterial colonies growing on solid medium.

ONPG broth (alternatively, commercially available prepared ONPG discs may be used according to the manufacturer's instructions).

Bacteriological straight wire/loop (preferably nichrome) or disposable alternative.

## 8 Quality control organisms

For Enterobacteriaceae,

**Positive control:**

*Escherichia coli*                      NCTC 10418 or NCTC 12241

**Negative control:**

*Proteus mirabilis*                      NCTC 10975

For *Neisseria* species,

**Positive control:**

*Neisseria lactamica*                      NCTC 10617

**Negative control:**

*Neisseria gonorrhoeae*                      NCTC 8375

**Note:** These strains are validated by NCTC to give this result.

## 9 Procedure and results

- a loopful of test organism from a culture plate or slant should be sufficient. Include the positive and negative controls with every batch of tests
- inoculate tubes containing ONPG reagent and incubate at 35-37°C for up to 24hr
- examine for yellow colour after 4hr and for up to 24hr

**Positive result:**


Yellow colour (indicates lactose fermenter).

**Negative result:**

Colourless/pale yellow (indicates lactose non-fermenter).



## Algorithm: ONPG ( $\beta$ -Galactosidase) test



**Note:**

**For Enterobacteriaceae**

**Positive control:** *Escherichia coli* NCTC 10418 or NCTC 12241

**Negative control:** *Proteus mirabilis* NCTC 10975

**For *Neisseria* species**

**Positive control:** *Neisseria lactamica* NCTC 10617

**Negative control:** *Neisseria gonorrhoeae* NCTC 8375

## References

An explanation of the reference assessment used is available in the [scientific information section on the UK SMI website](#).

1. Boadi S, Wren MW, Morris-Jones S. Selective testing of ss-galactosidase activity in the laboratory identification of Salmonella and Shigella species. J Clin Pathol 2010;63:1101-4. **B, III**
2. MacFaddin JF.  $\beta$ -Galactosidase (ONPG and PNPG) Tests. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 160-9. **B, III**
3. Lapage SP, Efstratiou A, Hill LR. The ortho-nitrophenol (ONPG) test and acid from lactose in Gram-negative genera. J Clin Pathol 1973;26:821-5. **B, II**
4. Clinical Microbiology Procedures Handbook: American Society for Microbiology; 2004. 3.3.2-3.3.2.13. **B, III**
5. Lapage SP, Jayaraman MS. Beta-Galactosidase and Lactose Fermentation in the Identification of Enterobacteria including Salmonellae. J Clin Pathol 1964;17:117-21. **B, II**
6. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2017-2018. 2017. **A, VI**
7. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices 1998. 1-37. **A, VI**
8. Home Office. Anti-terrorism, Crime and Security Act. 2001. **A, VI**
9. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books 2003. **A, VI**
10. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002 (as amended). Approved Code of Practice and guidance L5 (sixth edition). HSE Books,. 2013. **A, VI**
11. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 2009. **A, VI**
12. Health and Safety Executive. A Guide to Risk Assessment Requirements: Common Provisions in Health and Safety Law. HSE Books,. 2002. **A, VI**
13. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books,. 2002. **A, VI**

14. European Parliament. UK Standards for Microbiology Investigations (UK SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU in vitro Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes". 1998. **A, VI**
15. Department of Health. Transport of Infectious Substances. Best Practice Guidance for Microbiology Laboratories. Department of Health. 1-13. 2007. **A, VI**
16. Department for Transport. Transport of Infectious Substances, 2011 Revision 5. 2011. **A, VI**
17. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. MMWR Surveill Summ 2012;61:1-102. **B, V**
18. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 2005. 1-14. **A, VI**
19. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets 2000. **A, VI**
20. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive 2013. 1-35. **A, VI**
21. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive 2008. **A, VI**
22. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive 2005. **A, VI**
23. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office 2003. **A, VI**